

neaSNOM Microscope SDK




Overview

1. General design SDK / DLL
2. Setup access to automatically updated dll of neaSNOM
3. LabVIEW
 1. Connect to neaSNOM SDK with LabVIEW
 2. LabVIEW example: PsHet scan routine
1. Python
 1. Installation Python and pythonnet
 2. Python example: AFM scan routine
2. Documentation

1 General design neaSNOM SDK / DLL

The general design is depicted beneath, how to access any functionality of the neaSNOM:

1. **Connect to SDK**
2. **Create Microscope with the function ,Connect()' of the SDK**
=> You have now access to all functionallities of the AFM e.g. AutoApproach(), GotoTipPosition()
3. **You can now also connect to a Stage, Channel, or Scan:**
 - **Stage:** Allows you to move any stage in the Micscope, e.g. Parabolic mirror, sample
 - **Channel:** Allows you read out the current Value of any channel with CurrentValue
 - **Scan:** Allows you to create any Scan already implemented in the neaSNOM and you can set the parameter for the scan, e.g. AfmScan, WhitelightScan
4. **The function Start():**
 1. Starts the scan
 2. Creates one or all images of the channels available in this scan.



All SDK functionalities are documented under the following location: \\nea-server\updates\SDK\doc\index.html



2 Setup access to automatically updated dll of **neaSNOM**

Map Network Drive:

- 1) Click with the right mouse button on Computer in the Windows Explorer
- 2) Select Map network drive...
- 3) Name the Drive N: (any letter possible)
- 4) Connect to the following folder:

[\\nea-server\updates\SDK](#)


- 5) The SDK (DLL) in this loaction will be automatically updated with each **neaSCAN** update




3.1 Connect to neaSNOM with LabVIEW

neaspec LabVIEW SDK guide:

- 1) Start your LabVIEW Application
- 2) Open a new project or VI
- 3) Open the *Block Diagram View* by pressing 'Ctrl + E'
- 4) Right click on the white background and select:
Connectivity → *.NET* → *Constructor Node*
- 5) Place the constructor node in the block diagram
- 6) Select *Nea.Client.SDK* → *Connection*
- 7) To establish a connection type in the hostname via a *String Constant* and call it *nea-server*
- 8) Now open the .NET element *Invoke Node*
- 9) Wire it with the *Connection Node* and select *Connect()* as Method
- 10) The VI should look like this for now and the connection to the **neaSNOM** is established





2




3

8




6




nea-server

7



nea-server



9

3.2 LabVIEW Example: PsHet Scan Routine


After setting up the connection to the microscope all acquired Methods for standard Scan Routines and DAQ of the **neaSNOM** are accessible

Example PsHet Scan Routine:

- 1) Use another *Invoke Node* and wire it directly to *Connect* in the *Connection Node*
- 2) Choose the desired *Method* in the *Microscope Node*
- 3) In this example *PreparePsHetScan()* is selected as Method
- 4) Wire an *Invoke Node* to *PreparePsHetScan()*
- 5) This gives access to various scan parameters:

- Set Center
- Set Resolution
- Set Sampling Time
- Set Scan Area
- ... (further Settings in Documentation)

- 6) Connect all the parameters necessary for the scan in series
- 7) Wire the last *PsHetScan* Node with another Invoke Node and select the *Start* method to initiate the scan
- 8) Invoke the *Images 2D* Node to get access to individual channels in the *GetChannel* method
- 9) The *String Input* wired to *acronym* gives access to an user set Data Channel
- 10) For more than one channel, wire *GetChannel* methods in parallel
- 11) The path splits in two branches, data handling and scan status




3.2 LabVIEW Example: PsHet Scan Routine

Data Handling:

- 1) Within the *While Loop* the data flow for each pixel is saved to a customizable file location
- 2) In parallel the Image (*Image2D*) of the selected *Data Channel* is displayed in the *Front Panel*
- 3) All further data post processing may be done within the LabVIEW environment

Scan status indicators:

- 4) Display Scan Progress
- 5) LED to show, if the scan is completed



4.1 Installation Python and pythonnet

1. It is possible to use Python 2.7 and 3.7, but we recommend 3.7
2. Recommended download locations:
 - Anaconda includes already most important packages: <https://www.anaconda.com/distribution/#download-section>
 - Basic python (only basic packages installed): <https://www.python.org/downloads/>
3. After installations of python, install the packacke pythonnet:
“pip install pythonnet”

pythonnet allows python to handle the C# dll of the neaSNOM
4. Recommended packages to install (used in examples):
 - numpy
 - matplotlib
 - scipy
 - statistics

4.2 Python example: AFM scan routine

```
import sys
import time
import clr

#Basic steps to connect to Microscope
assembly_folder = '//nea-server/updates/SDK'
sys.path.append(assembly_folder)
clr.AddReference('Nea.Client.Hardware')
import Nea.Client.Hardware.SDK as neaSDK
neaClient = neaSDK.Connection('nea-server')
neaMic = neaClient.Connect()
time.sleep(0.1)

#Basic infos of Client and Server and print them to the console
ClVersion = neaMic.ClientVersion
print("Client Version: " + ClVersion)
SeVersion = neaMic.ServerVersion
print("Server Version: " + SeVersion)

#example to start auto approach after being connected
neaMic.CancelCurrentProcedure()
neaMic.RegulatorOff()
if not neaMic.IsInContact:
    neaMic.AutoApproach(setpoint)

#go out of contact and disconnect from the microscope
neaMic.RegulatorOff()
neaClient.Disconnect()
plt.show()
```

```
# sys class for importing folder
# time class for short delay
# pythonnet for reading C# SDK

#Import all DLLs in the folder or 'N:/updates/SDK' to mapped drive
# Import all DLLs in the folder
# Load the main DLL
# Import the DLL as element neaSDK
# Open up connection to microscope called neaClient
# Define the Microscope neaMIC
# Short delay makes things working fine

of python
# get Client Version

# get Server Version

# Cancel any running procedure (not obligatory)
# Retract sample (not obligatory)
# Check if system in contact (not obligatory)
# Start auto approach with chosen setpoint


# Retract sample
# Disconnect from SDK
```

This file `neaSNOM_AFMscan.py` is an example for a python script how to perform, show and save an AFM Scan with the SDK containing the example above!

5 Documentation

Overview of the whole SDK functionality is documented under the following location:

\\nea-server\updates\SDK\doc\index.html



Nea.Client.Hardware

Test Bla und so.

Nea.Client.Hardware.SDK Namespace

Type	Description
Connection	Provides a connection to remote neaSNOM controller.
Microscope	Represents a remote microscope.
Stage	Represents a stage with motors on three axes.

Nea.Client.Hardware.SDK.Images Namespace

Type	Description
Image2D	Represents a 2D image.
Image3D	Represents a 3D image.
Images<TScan,TImage>	Represents a collection of channels with scanned data.
Images2D	Represents a collection of 2D images.
Images3D	Represents a collection of 3D images.
Interferogram	Represents Interferograms.
Interferograms	Represents a collection of interferograms.
RawData	Represents a channel with scanned data.

Nea.Client.Hardware.SDK.Scan Namespace

Type	Description
AfmScan	Represents a scan with AFM parameters.
BasicScan	Represents a scan with basic parameters.
CfmScan	Represents a scan with CFM parameters.
FourierScan	Represents a scan with Fourier parameters.
PsHetScan	Represents a scan with PsHet parameters.
Scan2D	Represents a scan with 2D parameters.
Scan3D	Represents a scan with 3D parameters.
SerpentableScan2D	Represents a scan with basic parameters and optional serpent mode.
WhitelightScan	Represents a scan with Whitelight parameters.

Neaspec GmbH